

Euronda Monoart®

The original!

EN



Euronda®
Monoart

Euronda®

Over 40 years of
protection in dentistry

Euronda, synonymous with outstanding protection in the dental field for over 40 years. With an Italian origin and an international vocation, it has 4 branches operative and sells products in over 100 Countries worldwide. Euronda is one of the benchmark companies for infection control in the dental sector. Production includes two specific product lines, plus Euronda Alle® sterile disposable items.

Euronda Monoart®, a wide range of disposable products, accessories and personal protective equipment that teams maximum product quality with vibrant colours.

Euronda Pro System, the Euronda Total Protection System. A complete line of products to optimise the whole sterilisation process and make it safer. All Euronda products guarantee maximum quality and protection standards for dentists, professionals and patients.



Index

Monoart® Kit and Mix

4	Kit 5 Products
6	Kit Floral
8	Mix Saliva Ejectors
8	Mix Face Masks PTC 3
9	Mix Bandana
9	Mix Tray Paper
10	Mix Plastic Cups
11	Mix Apron PG30

Monoart® Disposable Line

12	Towel Up!
14	Saliva Ejectors EM15
15	Adapter for Saliva Ejector
16	EM19 Evo Aspirator Tips
18	EM21 Evo Aspirator Tips
19	EM40 Aspirator Tip
19	Adapter for Aspirator Tip
20	Disposable Latex Gloves
21	Disposable Nitrile Gloves
22	Plastic Cups
24	Plastic Cups Floral
25	Plastic Cups Dispenser
26	Face Masks Protection 3
28	Face Masks Protection 3 Floral
29	Face Masks Protection 4
30	Face Masks Protection 4 Sensitive
34	Bandana
36	Kit Dentist Chair Cover
37	Headrest Cover
38	Tray Paper
40	Apron PG30
42	Apron Dispenser
42	Apron PG20
43	Towel Dispenser
43	Baby Towel
43	Metal Chain for Towel
44	Instant Ice
45	Spray Ice
45	Cotton Rolls
45	Cotton Rolls Dispenser
46	Rubber Dam
46	Mixing Bowls
47	Facial Tissue
47	Bur Blocks
48	Kit Infection Control
48	Cut Proof Gloves
49	Gown
49	Apron
50	Barrier Film

Eurosteril®

51	Sterilization rolls
52	Mix Eurosteril®
52	Gusseted rolls
53	Sterilization pouches
53	Self-adhesive sterilization pouches

Monoart® Protective Glasses

56	Spheric Glasses
56	Light Glasses
57	Light Orange Glasses
57	Ultra Light Glasses
58	Evolution Glasses
58	Evolution Orange Glasses
59	Stretch Glasses
59	Stretch Orange Glasses
60	Baby Orange Glasses
60	Total Protection Glasses
61	Cube Glasses
61	Cube Orange Glasses
62	Small Operator Visor
62	Operator Visor





Kit Monoart® white



Kit Monoart® yellow



Kit Monoart® orange



Kit Monoart® lime



Kit Monoart® green



Kit Monoart® light blue



Kit Monoart® pink



Kit Monoart® lilac



Kit Monoart® blue

Kit Monoart® 5 Products

Practical Kit containing 5 different disposable products for 500 patients:

- 500 Monoart® Towel Up!
- 500 Monoart® EM15 Saliva ejector
- 500 Monoart® Plastic cup
- 200 Monoart® Face mask Protection 3
- 160 Monoart® Apron PG30 (Roll of 80 aprons)



Colours:

Lime, green, blue, light blue, lilac, pink, orange, yellow, white



Packaging:

1 carton containing 1 Kit





Kit Monoart® floral lime



Kit Monoart® floral light blue



Kit Monoart® floral pink



Kit Monoart® floral lilac

Kit Monoart® Floral

Monoart® Kit in floral colours for assistants and patients, containing:

- 100 Monoart® Plastic cups floral
- 50 Monoart® Face mask Protection 3 Floral
- 1 Monoart® Bandana
- 1 Monoart® Ultra Light Glasses



Colours:

Floral lime, floral light blue, floral lilac, floral pink



Packaging:

1 carton containing 1 Kit

Mix Monoart[®] Saliva Ejectors



Practical package containing 12 saliva ejector bags in 6 colours, two bags per colour.



Colours:

Blue, lime, yellow, orange, pink, lilac



Packaging:

Bag of 100 pieces
Carton of 12 bags

Mix Monoart[®] Face Masks Protection 3



Practical package containing 6 boxes of face masks Protection 3, one box per colour.



Colours:

Orange, yellow, lime, light blue, lilac, pink



Packaging:

Box of 50 pieces
Carton of 6 boxes

Mix Monoart® Bandana



Practical package containing 10 bandanas, one per colour.



Colours:

Lime, surgical green, light blue, blue, lilac, pink, orange, yellow, white, black



Packaging:

Carton of 10 multicoloured pieces

Mix Monoart® Tray Paper



Practical package containing 6 boxes of Monoart® Tray Paper (250 pcs).



Colours:

Yellow, orange, pink, lilac, light blue, lime



Packaging:

Box of 250 pieces
Carton of 6 boxes

Mix Monoart[®] Plastic Cups 1000



Practical package containing 1000 cups in 6 different colours (200 lime, 200 lilac, 150 orange, 150 green, 150 light blue, 150 pink).



Colours:

Green, lime, light blue, lilac, pink, orange



Packaging:

Bag of 50 pieces
Carton of 20 bags

Mix Monoart[®] Plastic Cups 3000



Practical package containing 3000 cups in 6 different colours (500 cups per colour).



Colours:


Lime, yellow, orange, pink, lilac, blue



Packaging:

Bag of 100 pieces
Carton of 30 bags

Mix Monoart[®] Apron PG30



Mix of 6 rolls of patient aprons PG30, one per colour.



Dimensions:

Rolls of 80 aprons
610 x 530 mm



Colours:

Lime, yellow, orange, pink,
lilac, blue



Packaging:

Carton of 6 rolls





Practical

Organize your space and keep your practice tidy, thanks to the compact design and the practical open/close packaging.



Monoart® Towel Up!



Monoart® Towel Up! disposable towels are made of pure cellulose which is highly absorbent. It is embossed and bonded with a waterproof polyethylene film. Thanks to the innovative packaging with a resealable open/close system, the towels are handled one at a time and they stay protected from dust and contamination, for maximum hygiene and highly practical use.

 **CE Type Ref :** MONOART® SALVIETTA, TOWEL, SERVIETTE



Classification:

Medical Device Class IIa



Colours:

Lime, emerald green, green, light blue, blue, lilac, pink, burgundy, orange, yellow, white, grey, black, blue lagoon



Dimensions:

330 x 450 mm



Packaging:

Box of 50 pieces
Carton of 10 boxes



Colourful

All the colours of the Monoart® line in a pack that expresses the taste and excellence of Made in Italy.



Safe

Absorbent, protective glue-free towels with fabric in pure cellulose that comes from a sustainable supply chain. Bonded using heat and embossing to a polyethylene film free of pollutants.



Hygienic

100% automated from the production line to the dental practice, the innovative open/close system reduces the risk of contamination.



Practical

The copper-clad wire is integrated in the tube to keep it at the set angle and ensure it won't protrude when bent.



Safe

Made of non-toxic PVC without DOP and phthalates, Class IIa Medical Devices with no expiry date.



Comfortable

With rounded fixed or removable tips that won't irritate mucous membranes.

Monoart[®] Saliva Ejectors EM15



Disposable, flexible saliva ejectors made of non-toxic, transparent PVC. The integrated metal wire perfectly holds them at the required angle and never protrudes during the bending process. Fixed or removable rounded tips are ideal for work on irritated mucous membranes.

 **CE Type Ref:** MONOART[®] ASPIRASALIVA EM15, EM15 SALIVA EJECTOR, EM15 SPEICHELSAUGER



Classification:

Medical Device Class IIa



Dimensions:

Diameter: 6.5 mm
Length: 12.5 or 15 cm



Colours:

Blue, green, lime, yellow, transparent, orange, pink, burgundy, lilac, blue lagoon



Packaging:

Bag of 100 pieces
Carton of 50 bags

Monoart[®] Adapter for Saliva Ejector



Safe, practical and autoclavable adapter for saliva ejector. Can be sterilized in autoclave at 134°C.

 **CE Type Ref:** ADATTATORE PER SALIVA, ADAPTER POR SALIVA EJECTORS, ADAPTER FÜR SPEICHELSAUGER



Classification:

Medical Device Class IIa



Packaging:

Bag of 10 pieces



Eco-friendly

Monoart® Aspirator Tips can go through up to 100 cycles in an autoclave at up to 134°C.



Comfortable

With a smaller contact surface thanks to the smaller end “brims” and fewer curves, the tips offer greater patient comfort.




Effective

The smooth interiors of the tips allow liquids to flow more easily.

Monoart[®] EM19 Evo Aspirator Tips



Monoart[®] EM19 Evo universal tips for surgical aspirators are recommended for procedures on adults. Ergonomic, non-slip grips make work easier for staff, while the linear shape keeps them quiet and promotes better hygiene and cleaning. EM19 Evo tips give outstanding patient comfort thanks to their rounded shape. They can be sterilized in autoclaves at up to 134°C.

 **CE Type Ref:** MONOART[®] CANNULA EM19,
EM19 ASPIRATOR TIPS, EM19 ABSAUGKANÜLEN



Classification:

Medical Device Class IIa



Dimensions:

Connection diameter:
16 mm; Length: 114 mm



Colours:

Lime, green, blue, lilac,
pink, burgundy, orange,
yellow, grey



Packaging:

Box of 10 pieces
of the same colour



Safe

The grooves are more extensive and reach all of the way to the end, for better grip.



Ergonomic

The tips are shorter and bent towards the end, for better control of operations on patients.



Silent

The angle of approximately 16° in the tips makes aspiration quieter.

Monoart[®] EM21 Evo Aspirator Tips



Monoart[®] EM21 Evo universal tips for surgical aspirators are recommended for procedures on cavities that are difficult to access. Ergonomic, non-slip grips make work easier for staff, while the linear shape keeps them quiet and promotes better hygiene and cleaning. They can be sterilized in autoclaves at up to 134°C.

 **CE Type Ref:** MONOART[®] CANNULA EM21, M21 ASPIRATOR TIPS, EM21 ABSAUGKANÜLEN



Classification:

Medical Device Class IIa



Dimensions:

Connection diameter:
16 mm; Length: 129 mm



Colours:

Lime, green, blue, lilac, pink, burgundy, orange, yellow, grey



Packaging:

Box of 10 pieces
of the same colour



Multisize

EM 19 Evo

Connection diameter: 16 mm
Length: 114 mm

EM 21 Evo

Connection diameter: 16 mm
Length: 129 mm

EM 40

Connection diameter: 11 mm
Length: 200 mm



Monoart® EM40 Aspirator Tip



Universal moulded tip for surgical aspirator.
Can be sterilised in autoclave at 134°C.

 **CE Type Ref:** MONOART® CANNULA EM40,
EM40 ASPIRATOR TIPS, EM40 ABSAUGKANÜLEN



Classification:

Medical Device Class IIa



Dimensions:

Connection diameter:
11 mm; Length: 200 mm



Colour:

Grey



Packaging:

Box of 10 pieces

Monoart® Adapter for Aspirator Tip



Safe and practical adapter for EM40 aspirator tip.
Can be sterilised in autoclave at 134°C.

 **CE Type Ref:** MONOART® ADATTATORE PER CANNULA,
ADAPTER POR ASPIRATOR TIPS, ADAPTER FÜR ABSAUGKANÜLEN



Classification:

Medical Device Class IIa



Packaging:

Bag of 10 pieces

Monoart[®] Disposable Latex Gloves



Safe, comfortable and durable: Monoart[®] disposable latex gloves provide hygienic, comfortable protection for hands. They are easy to wear, stretchy and guarantee outstanding tactile sensitivity. In addition, the rough micro-texture on the surface gives superb grip in all conditions, so they are ideal for working with complete peace of mind during operations of all kinds. Monoart[®] latex gloves come in a range of bright colours to give a unique look to every dental practice.

 **CE Type Ref :** MONOART[®] LATEX GLOVES



Classification:

Medical Devices Class I, pursuant to EN 455-1, EN 455-2, EN 455-3 and EN 455-4



Weight:

6.2 ± 0.3 gr



Available sizes:

XS, S, M, L



Colours:

Lime, blue, lilac, pink, white



Packaging:

Box of 100 gloves
Carton of 10 boxes



Safe

They provide an exceptional protective barrier against contamination and infection and they do not contain lubricant powders, chemical residue and endotoxins, which are often responsible for the respiratory symptoms caused by latex.



Certified

They comply with strict quality control criteria:

- AQL 1.5
- Level G1
- Low levels of protein residue (< 50 mg/g).

They are certified in accordance with the international directive 93/42/EEC, as amended by directive 2007/47/EC.



Durable

They are water, acid and alkali resistant and they provide good protection against micro-organisms and the majority of caustic substances and cleaning products.

Monoart[®] Nitrile Gloves



The Monoart[®] disposable nitrile gloves, conceived as an alternative for users allergic to latex proteins, provide an exceptional protective barrier against contamination and infections. They offer extraordinary mechanical and chemical resistance, making them suitable for all those situations in which there is a risk of perforation, cuts, piercing or handling particularly aggressive chemicals. The quality materials also guarantee excellent tactile sensitivity, while the micro-rough texture ensures excellent grip in all conditions.

 **CE Type Ref:** MONOART[®] NITRILE GLOVES



Classification:

Medical Devices Class I,
pursuant to EN 455-1,
EN 455-2, EN 455-3
and EN 455-4



Available sizes:

XS, S, M, L



Colours:

Black



Packaging:

Box of 100 gloves
Carton of 10 boxes





Ergonomic

With a rounded edge and a secure grip.



Capacious

Available in 166 ml (white only) and 200 ml (multicolour) versions.



Resistant

Made from a certified, hygienic and practical plastic material.

Monoart[®] Plastic Cups Coloured



Disposable plastic cups with rounded edges and a capacity of 200 cc made from non-toxic polypropylene. Monoart[®] disposable cups coloured have been made with a new material and the polystyrene has been replaced with polypropylene. All the new polypropylene cups will bear an embossed company logo, testifying to their high quality and safety standards Euronda. Available in numerous colours to customize your dental practice.



Capacity:
200 cc



Colours:
Light blue, blue, lilac, pink, fuchsia, red, orange, yellow, lime, green, black, blue lagoon



Packaging:
Bag of 100 pieces
Carton of 30 bags

Monoart[®] Plastic Cups White



Disposable white plastic cups made from non-toxic polypropylene with rounded edges and a capacity of 166 cc.



Capacity:
166 cc



Colour:
White



Packaging:
Bag of 100 pieces
Carton of 30 bags



Monoart[®] Plastic Cups Floral



200 cc disposable plastic cups with rounded edges and a floral pattern. Bring a blast of colour, design and vibrancy to your dental practice with their floral pattern. Disposable floral cups are made with a new material: polypropylene instead of polystyrene. Unaltered capacity and ergonomics, but enhanced resistance and reliability.



Capacity:
200 cc



Colours:
Floral lime, floral light blue,
floral pink, floral lilac



Packaging:
Bag of 100 pieces
Carton of 30 bags

Monoart[®] Plastic Cups Dispenser



Wall-mounted handy dispenser for Monoart[®] plastic cups. It contains up to 70 coloured Monoart[®] plastic cups with a capacity of 200 ml, or up to 90 white Monoart[®] plastic cups with a capacity of 166 ml.



Colour:
White with transparent
smoky shade body



Safe

Available in two versions: with 3 or 4 layers. They are made entirely of hypo-allergenic materials and guaranteed to last for 3 hours.



Comfortable


Absorbent non-woven fabric lining for effective absorption of perspiration and humidity. Built-in nose bridge.



Protective

External layer made of water-repellent non-woven fabric for protection from dust and moisture. Lining with particle and bacteria filters. Protection 4 also comes with a polymeric filter for excellent resistance against fluids.

Monoart[®] Face Masks Protection 3



Disposable, hypo-allergenic face mask with 3 layers:

- Soft, coloured and water-repellent non-woven fabric for improved protection against dust and humidity;
- PFE (Particle Filtration Efficiency) filter and BFE (Bacterial Filtration Efficiency) filter > 99%;
- Comfortable non-woven cotton lining for extra absorbance of sweat and humidity.

It is ideal for protection from dust and aerosols: $\Delta P < 3.0 \text{ mm H}_2\text{O/cm}^2$ in accordance with the EN 14683 standard (differential pressure: surgical mask breathability rating). Built-in nose bridge.



Classification:

Medical Device Class I
EN 14683: Type II
ASTM F2100: Level I



Duration:

3 hours of protection



Colours:

White, yellow, lime, green, light blue, lilac, pink, orange



Packaging:

Box of 50 masks with an "open/close" system.
Carton of 20 boxes



CE Type Ref: Monoart[®] FACE MASK PROTECTION 3



Soft coloured and water repellent non-woven fabric

Particle and bacterial filters
PFE and BFE > 99%

Comfortable non-woven cotton

Monoart[®] Face Masks Protection 3 Floral



Disposable hypo-allergenic face mask, composed by 4 layers:

- Soft, coloured and water-repellent non-woven fabric with floral motifs;
- PFE (Particle Filtration Efficiency) filter and BFE (Bacterial Filtration Efficiency) > 99%;
- 2 layers of comfortable, absorbent non-woven fabric lining for effective absorption of sweat and moisture.

It is ideal for protection from dust and aerosols:

$\Delta P < 3.0 \text{ mm H}_2\text{O}/\text{cm}^2$ in accordance with the EN 14683 standard (differential pressure: surgical mask breathability rating). Built-in nose bridge.



Classification:

Medical Device Class I
EN 14683: Type II
ASTM F2100: Level I



Duration:

3 hours of protection



Floral colours:

Floral lime, floral light blue, floral lilac, floral pink, floral burgundy




Packaging:

Box of 50 masks with an "open/close" system.
Carton of 20 boxes



CE Type Ref: Monoart[®] FACE MASK PROTECTION 3



Monoart[®] Face Masks Protection 4



Disposable, hypo-allergenic face mask with 4 layers:

- Soft, colourful water-repellent non-woven fabric for unparalleled protection from dust and moisture;
- PFE (Particle Filtration Efficiency) filter and BFE (Bacterial Filtration Efficiency) filter > 99%;
- Polymeric filter for excellent fluid resistance, 120 mmHg;
- comfortable non-woven cotton lining for extra absorbance of sweat and humidity when used by operators.

It is ideal for protection from dust and aerosols: $\Delta P < 3.0 \text{ mm H}_2\text{O}/\text{cm}^2$ in accordance with the EN 14683 standard (differential pressure: surgical mask breathability rating). Built-in nose bridge.



Classification:

Medical Device Class I
EN 14683: Type II
ASTM F2100: Level I



Duration:

3 hours of protection



Colours:


Green, light blue,
orange, black



Packaging:

Box of 50 masks with an
"open/close" system.
Carton of 20 boxes

 **CE Type Ref:** Monoart[®] FACE MASK PROTECTION 4





Comfortable

Internal lining in hydrophilic NWF, Perfect-Fit system, externally welded elasticated ties and an elongated nose bridge for an extremely close comfortable fit.



Protective

External layer in water-repellent non-woven fabric against dust and moisture. Internal layer with particulate and bacterial filter and hydrophobic filter for outstanding fluid resistance.



Latex Free

Devoid of latex, glass fibre, pollutants and chemical substances, they have been made from hypo-allergenic cellulose for sensitive skins.

Monoart[®] Face Masks Protection 4 Sensitive



The Monoart[®] PTC4 Sensitive Face Masks are class I disposable medical devices and have been developed to offer total operational protection (the same standards as the Protection 4 range), but they also deliver the safety levels of hypo-allergenic materials which have been carefully sourced to respect sensitive skins and avoid the risk of irritation and dermatitis.

Resistant face mask, ideal for protection against dust, aerosol, humidity and splashes:

- Soft and colourful water-repellent NWF for optimal protection against dust, powder and residue.
- A particulate filter (PFE, Particle Filtration Efficiency) and a bacterial filter (BFE, Bacterial Filtration Efficiency) > 99%.
- Improved hydrophobic layer for excellent fluid resistance: 160 mm Hg.
- Comfortable hydrophilic non-woven fabric inside which effectively absorbs sweat and moisture.



Classification:

Medical Device Class I
EN 14683: Type II
ASTM F2100: Level IIR



Color:

Light Blue




Packaging:

Box of 50 masks with an
“open/close” system.
Carton of 20 boxes



CE Type Ref: Monoart[®] FACE MASK PROTECTION 4
SENSITIVE



Soft coloured and water repellent non-woven fabric

Polymeric fluid resistant filter

Particle and bacterial filters PFE and BFE > 99%

Comfortable hydrophilic non-woven fabric

Monoart[®] Face Mask certifications

All the Monoart[®] Face Masks are class I medical devices, compliant with the ASTM F 2100 standards, according to the EN 14683 standard. Monoart[®] Protection 3 and Protection 3 Floral Masks are Type II, while Monoart[®] Protection 4 and Protection 4 Sensitive Masks are of type IIR.

The resistance of the Monoart[®] Face Masks to the penetration of splashing liquids has been tested according to the ASTM F1862 standard, and is compliant, as requested by the European standard EN 14683, with the values below:

Test	Type I	Type II PTC3 and PTC3 Floral	Tipo IIR PTC4 and PTC4 Sensitive
Bacterial filtration efficiency (BFE) (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 29,4	< 29,4	< 49,0
Splash resistance pressure (mm Hg)	Not required	Not required	≥ 120

Why use the face mask?

The face mask offers protection from the chemical substances used during surgery and reduces the risk of contamination from viruses and bacteria that can cause infections in both professionals and patients.

How to use the face mask correctly?

The mask must adhere well to the face and remain dry at all times.

When to change the face mask

The mask must be changed if it is exposed to nebulised substances, or if it gets wet, and in any case for each patient treated.





Monoart[®] Bandana



Pleasant and comfortable to wear, this cotton bandana has adjustable ties and can fit all sizes. Thanks to the top-quality, colour-fast cotton used, it is autoclavable at 121°C. Choose the ideal colour to match your style.

 **CE Type Ref :** MONOART[®] BANDANA



Classification:

Personal Protective Equipment



Colours:

Lime, green, lilac, surgical green, light blue, blue, pink, orange, yellow, white, black



Packaging:

Single pack
Carton of 10 pieces
of the same colour



Comfortable

Cotton, comfortable to wear bandanas.



Autoclavable

Autoclavable at 121°C thanks to the top-quality cotton used.



Fits All Sizes

They can fit all sizes thanks to their adjustable ties.



Comfortable

Made of soft, hypo-allergenic and breathable non-woven polypropylene fabric.



Hygienic

Daily dentist chair protection for patient safety.



Universal

A headrest cover, backrest cover and seat cover are all included. It can be adjusted to fit all dentist chairs thanks to the stretchy edges and the ties.



Kit Monoart[®] Dentist Chair Cover



A disposable universal dentist chair cover for the protection and safety of patients. It is made of soft non-woven polypropylene fabric, which is hypo-allergenic and breathable. The Kit contains a headrest cover, a backrest cover and a seat cover that are all made of 30 g/m² non-woven fabric. Its universal size perfectly fits all dentist chair models thanks to the stretchy edges and the ties, which make it extremely versatile.



Colours:

White, surgical green, light blue, blue, lime, yellow, orange, lilac



Packaging:

Bag of 1 dentist chair cover Kit
Carton of 25 single colour bags

Monoart[®] Headrest Cover



Polyethylene headrest cover, with cellulose exterior and polyethylene lining. In a practical dispenser.



Colours:

White, light blue, green



Dimensions:

280 x 300 mm



Packaging:

Box of 250 pieces
Carton of 8 boxes



Monoart[®] Tray Paper



Chlorine-, dioxin- and ash-free crinkled paper that is permeable to steam and gas. Thanks to its advanced technical features, steam does not disintegrate or weaken this paper, so it is ideal for lining trays.



Colours:

White, yellow, orange, pink, lilac, light blue, green, lime



Dimensions:

280 x 180 mm



Packaging:

Box of 250 pieces
Carton of 10 boxes



Practical

The ideal base for the bottom of a tray.



Resistant

Permeable to steam and gas.



Safe

Chlorine- and dioxin-free.



Ergonomic

The innovative system with laces, allows to tie the apron directly to the patient's neck.



Colourful

They come in numerous colours that can be matched with other Monoart® products.



Protective

Practical shape that perfectly covers the patient's neck and shoulders.



Monoart[®] Apron PG30



A roll of disposable aprons for patient protection. They are made of a layer of highly absorbent paper and a layer of waterproof polyethylene, which are bonded and embossed. Unlike traditional towels, they provide perfect protection for patients thanks to their large size and the specially shaped neck and shoulder sections.



Dimensions:

Roll of 80 aprons:
610x530 mm

Roll of 60 aprons:
810x530 mm



Colours:


Lime, green, emerald green, blue, light blue, lilac, pink, burgundy, orange, yellow, white, blue lagoon



Packaging:

Carton of 6 rolls

Instructions for use



Hygienic

No glue is used. Embossing and heating are used to bond the polyethylene film and fabric together.



Absorbent

It guarantees high speed of absorbency to protect the patient.

Monoart[®] Apron dispenser



Practical dispenser for aprons on a roll.




Colour:
White RAL 9002



Dimensions:
343 x 132 x 142 mm

Monoart[®] Apron PG20

Roll of patient aprons made from
a waterproof polyethylene layer.



Dimensions:
Rolls of 200 aprons
810 x 530 mm



Colours:
White, green, blue



Packaging:
Carton of 6 rolls

Monoart® Towel dispenser



Practical dispenser for disposable towels.



Colour:
White RAL 9002



Dimensions:
165 x 165 x 125 mm

Monoart® Baby towel



Pure padding of highly absorbent cellulose coupled with a polyethylene waterproof film with cheerful designs.



Dimensions:
350 x 400 mm



Colours:
Light blue, pink



Packaging:
Roll of 100 bibs
Carton of 24 rolls

Monoart® Metal chain for towel



Adjustable metal chain for towel.



Length:
Ca. 45 cm



Packaging:
Single pack

Monoart[®] Instant Ice



Thanks to its formula, Euronda Monoart[®] K2P-Cellulose is one of the best products of its kind on the market in terms of the cold produced, duration, functionality, patient safety and the aesthetic appeal of the packaging. The new mixture of nitrogenous salts guarantees better quality and longer lasting cooling, while the smaller size of the bag gives greater concentration in the treatment zone and takes up less space.

 **CE Type Ref:** MONOART[®] GHIACCIO ISTANTANEO,
INSTANT ICE KÄLTEKOMPRESSE



Classification:

Medical Device Class IIa



Size:

150 x 175 mm



Packaging:

Carton of 24 bags



Shelf life:

5 years



Comfortable

The new material in K2P-Cellulose guarantees a better tactile sensitivity and pleasant effect when on contact with the patient's skin.



Practical

The smaller size of the bag concentrates the effects on the relevant treatment zone.



Lasting

It is pressure activated, the new mixture of nitrogenous salts guarantees better quality and longer lasting cooling.

Monoart[®] Spray Ice



Cooling spray in a 200 ml aerosol. Ideal for cooling materials (impression, etc...)



Packaging:

Carton of 12 cans



Shelf life:

3 years

Monoart[®] Cotton Rolls



Long fibre, pure cotton rolls that are compact and highly absorbent. Cellulose-free. It is a Medical Device Class I.

Ø 8 mm / L 37 mm:

Box of 300 g

Carton of 30 boxes

Ø 10 mm / L 37 mm:

Box of 300 g

Carton of 60 boxes

Ø 12 mm / L 37 mm:

Box of 300 g

Carton of 30 boxes



CE Type Ref: MONOART[®] RULLI DI COTONE, COTTON ROLLS, WATTERROLLEN

Monoart[®] Cotton Rolls Dispenser



A practical dispenser for cotton rolls.



Colours:

White, light blue

Monoart[®] Rubber Dam



Pure latex rubber dam, medium thickness.

 **CE Type Ref:** DIGA DENTALE, RUBBER DAM, KOFFERDAMM



Classification:

Medical Device Class I



Colours:

Grey, green



Packaging:

Box of 36 pieces

Mixing Bowls



A varied selection of green mixing bowls.

Mixing bowl n.1 (small)

Diameter: 110 mm
Height: 75 mm

Mixing bowl n.2 (medium)

Diameter: 125 mm
Height: 90 mm

Mixing bowl n.3 (large)

Diameter: 140 mm
Height: 112 mm

Mixing bowl n.4 (maxi)

Diameter: 165 mm
Height: 154 mm

Monoart[®] Facial Tissue



Facial tissue made from two plies of pure cellulose paper, soft and resistant.



Dimensions:

200 x 210 mm



Colour:

White



Packaging:

Single pack (100 tissues)
Carton of 40 boxes

Bur Blocks



Bur blocks for dentists with autoclavable base.

Bur blocks with autoclavable base for:
18FG, 18RA or
18 FG RA mixed.

Bur blocks with autoclavable base for:
36FG, 36RA or
36 FG RA mixed.

Bur blocks with autoclavable base for:
72FG, 72RA or
72 FG RA mixed.

Kit Monoart® Infection Control



A safe barrier against infections.

Practical Kit containing:

- 50 Monoart® Apron
- 50 Monoart® Face mask Protection 3
- 2 Monoart® Bandana
- 2 Monoart® Ultra Light glasses
- 2 Monoart® Gowns
- 1 pair of cut proof gloves



Packaging:

1 carton containing 1 Kit

Cut proof gloves



Antiallergic gloves with anti-slip surface and a shaped form for higher comfort, sensitivity and protection for the user. It stands punctures, and protects against aggressive chemicals, solvents and oils. With linen nitrile, it can be sterilized in the autoclave and it is latex free.



Available sizes:

S, M



Colour:

Violet



Packaging:

Bag of 1 pair
Carton of 9 bags

Monoart[®] Gown



Gown with zip, mandarin collar and elastic cuffs for the total protection of the operator.



Classification:

Personal Protection Equipment



Colour:

White



Packaging:

Carton of 25 pieces

Monoart[®] Apron



Practical disposable polyethylene apron with a comfortable fit to protect operators.



Colour:

White



Dimensions:

80x110 cm



Packaging:

Bag of 50 pieces
Carton of 9 bags

Monoart[®] Barrier film



Disposable film for the protection of areas and surfaces, especially those difficult to reach, clean or disinfect. The light adhesive of the barrier film adheres easily to surfaces without leaving any residue. The barrier film has a non-adhesive edge on both sides, to simplify placement and removal. The use of the barrier film reduces the risk of cross-contamination caused by aerosol and spraying. The pack contains 1 roll of 1200 tear-off 100x150 mm sheets with dispenser.

**Classification:**

DM Class I

**Dimension:**

100 x 150 mm

**Color:**

Blu

**Packaging:**

Single pack

 **CE Type Ref:** MONOART[®] BARRIER FILM



Eurosteril® Sterilization roll



Eurosteril® sterilization rolls are made of the best certified materials: heavy weight medical paper (60 g/m²) and a double layer of light blue polyester/polypropylene film. This gives them a guaranteed seal during sterilization, makes them extremely easy to open without releasing fibres or tearing, and offers maximum microbial protection.

If stored properly, they can be sterile for up to 6 months. In addition, the chemical colour indicators signal a successfully completed cycle in the autoclave, thus guaranteeing maximum safety for patients.

 **CE Type Ref:** EUROSTERIL® ROTOLI PER STERILIZZAZIONE, STERILIZATION ROLLS, STERILISATIONS ROLLEN



Classificazione:

DM Class I products fully compliant with Directive 93/42/EEC (as amended by 2007/47/CE) and the EN 868-5 standard.



Packaging:

5 cm x 200 m
Carton of 8 rolls

7.5 cm x 200 m
Carton of 5 rolls

10 cm x 200 m
Carton of 4 rolls

15 cm x 200 m
Carton of 2 rolls

20 cm x 200 m
Carton of 2 rolls

25 cm x 200 m
Carton of 1 roll

30 cm x 200 m
Carton of 1 roll

40 cm x 200 m
Carton of 1 roll



Long-lasting

They offer a guaranteed seal during sterilization and keep items sterile for up to 6 months in suitable storage conditions and subject to validation of the sterilisation and sealing process.



Safe

Excellent peelability during opening, with no release of paper fibres or tearing in the plastic material.



User-friendly

The Class I chemical indicators change colour during the sterilization process.

Mix Eurosteril®



Practical Mix of rolls made of medical high weight paper (60g/m²) and double layer film of polyester/polypropylene, maintaining the durability of the instrument's sterility up to 6 months.



Mix composed by:

2 rolls 5 cm
2 rolls 7.5 cm
1 roll 10 cm



Packaging:

Carton of 1 mix

Eurosteril® Gusseted rolls



Gusseted rolls made of medical high weight paper (60g/m²) and double layer film of polyester/polypropylene, maintaining the durability of the instrument's sterility up to 6 months.



Classification:

Medical Device Class I



Packaging:


20 cm x 5 x 100 m
Carton of 2 rolls

25 cm x 5 x 100 m
Carton of 1 roll

25 cm x 6,5 x 100 m
Carton of 1 roll

 **CE Type Ref:** EUROSTERIL® ROTOLI PER STERILIZZAZIONE,
STERILIZATION ROLLS, STERILISATIONS ROLLEN

Eurosteril® Sterilization pouches



Thermal-welding standard sterilization pouches made of heavy weight (60 g/m²) white medical paper coupled with a light blue polyester/polypropylene layer. Three sided weldings with impermeable and uniform canals for the highest microbe protection. The changing chemical colour indicators verify the sterilization process by steam or EtO gas (ethylene oxide). Appropriate peeling without the release of paper fibres or laceration of the film at the moment of opening the packaging.



Classification:

Medical Device Class I



Packaging:

5 x 25 cm
Box of 500 pieces
7,5 x 25 cm
Box of 500 pieces
10 x 25 cm
Box of 500 pieces



CE Type Ref: EUROSTERIL® BUSTE PER STERILIZZAZIONE, STERILIZATION POUCHES, STERILISATIONS BEUTEL

Eurosteril® Self-adhesive sterilization pouches



Self adhesive sterilization pouches made of heavy weight (60 g/m²) white medical paper coupled with a light blue polyester/polypropylene layer. Three sided weldings with impermeable and uniform canals for the highest microbe protection. Adhesive band closure of synthetic rubber. The changing chemical colour indicators verify the sterilization process by steam or EtO gas (ethylene oxide). Appropriate peeling without the release of paper fibres or laceration of the film at the moment of opening the packaging.



Classification:

Medical Device Class I



Packaging:

9 x 25 cm
Box of 200 pieces
Carton of 9 boxes
14 x 26 cm
Box of 200 pieces
Carton of 6 boxes
19 x 33 cm
Box of 200 pieces
Carton of 5 boxes



CE Type Ref: EUROSTERIL® BUSTE PER STERILIZZAZIONE AUTOSIGILLANTI, SELF ADHESIVE STERILIZATION POUCHES, STERILISATIONSBEUTEL SELBSTKLEBEND



Comfortable

The patented SoftPad frame stems (for the Cube, Cube Orange and Stretch models) offer greater protection, stability and adaptability to the face and can be adjusted to three different lengths. The lenses in the orange models reduce eye fatigue.



Resistant

Permanent Univet Double Coating prevents fogging inside the lenses (pursuant to requirement N of the EN166 standard) and scratching outside (pursuant to requirement K of the EN 166 standard). No more scratches or fogging, for outstanding performance.



Ergonomic

They are extremely light and have non-slip nose bridges. Versions for fitting over corrective glasses are also available.



Monoart[®] Protective Glasses

Monoart[®] protective glasses for staff and patients are designed and made to combine maximum safety in every situation with exceptional comfort and stylish designs.



Anti Scratch

Basic Anti-scratch treatment protects against scratches and abrasions.

*K Anti-scratch treatment



Anti Scratch Plus

Anti-Scratch Plus treatment constitutes a barrier against surface damage, it protects against scratches and abrasions. It is certified according to the requirement of EN166 "K".



Anti Fog

Basic Anti-fog treatment that prevents the formation of condensation.

*N Anti-fog treatment



Anti Fog Plus

Anti-Fog Plus treatment is a permanent treatment that prevents the formation of condensation. It is certified according to the requirement of EN166 "N".



Univet Double Coating

The revolutionary coating created through UDC technology provides to your eyewear superior resistance from scratches and fogging for outstanding performance in any working condition. Deposition of coatings in separate processes on the surfaces of the lens:

- exclusive anti-fog coating inside the lens
- multilayer scratch-resistant coating outside the lens. Both treatments are permanent and assure to the product an highly superior performance to EN 166 standards. This technology is expressed in Monoart[®] Stretch Orange glasses.



SoftPad Technology

SoftPad is a patented terminal that has been added to the legs of the eyewear to offer comfort and protection. Thanks to the innovative use of materials and a meticulous design, the ergonomics of these glasses reach a new level of comfort, providing perfect stability to the safety glasses and full adaptability to the face of the user. In addition to the length adjustment in three positions, the SoftPad system has a further tilting movement that allows the device to perfectly fit the user's head. The SoftPad system is supplied with the Monoart[®] Cube glasses, Cube Orange glasses and Stretch glasses.

Monoart[®] Spheric Glasses



Glasses for staff and patients with a single, optically controlled spherical polycarbonate lens that offers a 180° field of view. The length and angle of the temples are adjustable. They are made from over-moulded material that guarantees excellent fit and comfort. Polycarbonate nose bridge. Neutral, colourless anti-scratch and anti-fog lenses protect against impacts and UV rays.



Classification:
Personal Protection
Equipment



Colour:
White and black frame



Weight:
36 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Light Glasses



Protective glasses for patients or operators, basic line with air ventilation system incorporated in the side protections. Single lens in optically checked polycarbonate: these protective glasses can be fitted over corrective glasses. Their transparent colourless lens protects against impact and UV rays to guarantee optimal visibility conditions.



Classification:
Personal Protection
Equipment



Colour:
Transparent frame



Weight:
43 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Light Orange Glasses



Over-glasses for patients and staff with a single polycarbonate orange lens which completely covers corrective glasses. The length and angle of the temples can be adjusted. The temples are made from over-moulded material and they guarantee excellent fit and comfort. The orange lens with anti-fog and anti-scratch treatment protects against impacts, UV rays and blue light (525 nm). It improves contrast and enhances visual sharpness to reduce eye fatigue. Ideal for use during dental composite polymerization.



Classification:
Personal Protection
Equipment



Colour:
Black and orange frame



Weight:
43 g



Packaging:
Single piece
Carton of 10 pieces



Monoart[®] Ultra Light Glasses



Protective glasses for staff and patients, with ergonomic, adjustable length temples, protection incorporated within the lens, a polycarbonate nose bridge and eyebrow protection. The colourless lens protects against impacts and UV rays to guarantee optimum visibility.



Classification:
Personal Protection
Equipment



Colour:
Pink frame



Weight:
29 g



Packaging:
Single piece
Carton of 10 pieces



Monoart[®] Evolution Glasses

Glasses for staff and patients with an extremely lightweight, modern design. The temples guarantee excellent fit and comfort. They are made from over-moulded material and feature a memory system that keeps them in the desired shape. The wrap-around transparent lens is made of heavy-duty polycarbonate that protects against impacts and UV rays. They are equipped with a non-slip rubber nose bridge for maximum comfort.

**Classification:**

Personal Protection Equipment

**Colour:**

White and black frame

**Weight:**

25 g

**Packaging:**

Single piece
Carton of 10 pieces



Monoart[®] Evolution Orange Glasses

Glasses for staff and patients with a lightweight, modern design. The temples guarantee excellent fit and comfort. They are made from over-moulded material and feature a memory system that keeps them in the desired shape. The wrap-around orange lens is made of heavy-duty polycarbonate that protects against impacts, UV rays and blue light (525 nm). It improves contrast and enhances visual sharpness while reducing eye fatigue. Ideal for use during dental composite polymerization. Non-slip rubber nose bridge. Can fit over corrective glasses.

**Classification:**

Personal Protection Equipment

**Colour:**

Black and orange frame

**Weight:**

25 g

**Packaging:**

Single piece
Carton of 10 pieces



Monoart® Stretch Glasses

Stretch glasses with a frame designed to buffer the contact points between the glasses and the face, while also making them lighter and more unobtrusive for users. The sporty, wrap-around design features adjustable length temples that are ultra-flexible and extremely twist proof, for an unbeatably comfortable fit. The heavy-duty, optically controlled transparent lens with anti-fog and anti-scratch treatment protects against UV400 rays (sunlight).



Classification:
Personal Protection
Equipment



Colour:
White and blue frame



Weight:
26 g



Packaging:
Single piece
Carton of 10 pieces

Monoart® Stretch Orange Glasses

Orange glasses with a frame designed to buffer the contact points between the glasses and the face. The ultra-flexible, sporty design is extremely twist proof and gives an unbeatably comfortable fit. There is an eyebrow ventilation system and a wrap-around lens to give extra protection around the eyes. The heavy-duty, optically controlled orange lenses with anti-fog and anti-scratch treatment protect against UV400 rays and UV525 rays (polymerization lamps). A soft rubber nose bridge keeps the glasses stable.



Classification:
Personal Protection
Equipment



Colour:
Black and orange frame



Weight:
31 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Baby Orange Glasses

Lightweight, easy-to-manage glasses for children. 100% polycarbonate, with black temples and an orange lens that protects against UV and UV525 (blue light) rays. Basic anti-fog treatment.



Classification:
Personal Protection
Equipment



Colour:
Black and orange frame



Weight:
23 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Total Protection Glasses

Glasses for staff and patients with an optically controlled lens and built-in side protection. The length and angle of the temples are adjustable. They are made from over-moulded material to guarantee excellent fit and comfort. Rubber eyebrow protection and non-slip soft rubber nose bridge. Neutral, colourless anti-scratch and anti-fog lenses protect against UV rays and guarantee both optimum visibility and excellent impact resistance.



Classification:
Personal Protection
Equipment



Colour:
Grey and coral frame



Weight:
45 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Cube Glasses



Transparent operator over-glasses, with an extremely modern and wrap-around design. Adjustable frame stems guarantee an excellent fit and comfort, and memorise the desired shape of the user thanks to the SoftPad technology. Transparent optically checked polycarbonate lens, providing a 180 degree view. Heavy duty lens, with anti-fog and anti-scratch treatment, that protect against impact and UV400 (sunlight) rays. Can be positioned over any corrective glasses.



Classification:
Personal Protection
Equipment



Colour:
White and blue frame



Weight:
37 g



Packaging:
Single piece
Carton of 10 pieces

Monoart[®] Cube Orange Glasses



Orange over-glasses for staff with a cutting-edge design featuring a wrap-around square shape. Adjustable length temples guarantee an excellent fit and comfort, while the SoftPad technology keeps them in the user's desired shape. Orange optically checked polycarbonate lens, providing a 180° view. The heavy-duty lens with anti-fog and anti-scratch treatment protects against impacts, UV400 (sunlight) rays and UV25 blue light rays (polymerization lamps). They can fit over any corrective glasses.



Classification:
Personal Protection
Equipment



Colour:
Black and orange frame



Weight:
37 g



Packaging:
Single piece
Carton of 10 pieces

Small Operator Visor



Disposable short operator visor designed to cover the eyes and cheekbones. Ultra-lightweight. Replaceable transparent neutral shields. Frames included in a pack: yellow, blue, green, pink and black.



Classification:
Personal Protection Equipment



Colours:
Pink, yellow, green, blue and black frames



Weight:
9 g



Packaging:
Carton of 5 visors with 20 shields

Monoart[®] Operator Visor



Disposable operator visor with transparent frame and replaceable shield. Transparent and anti-fog disposable shield. Ultra light visor.



Classification:
Personal Protection Equipment



Colour:
Transparent frame



Weight:
25 g



Packaging:
Box of 10 visor with 10x10 shields or refill with 20 shields




Eyewear marking

Added to personal protective equipments (PPE) the marking guarantees respect of current safety regulations. Required by the law, the marking is added both on the lens, to identify scale number, manufacturer, optical class, mechanical resistance and optional requirements, and on the frame to identify manufacturer, reference normative and mechanical resistance.

European Normatives

EN 166 Personal eye-protection - Specifications
 EN 167 Personal eye-protection - Optical test methods
 EN 170 Personal eye-protection - Ultraviolet filters.
 Transmittance requirements and recommended use.

Euronda is the exclusive distributor of Univet products in the dental field.



Frame marking identification

U	EN 166	FT	CE
Manufacturer	Standard	Mechanical resistance	

Lens marking identification

2 - 3	U	1	FT	KN*	CE
Scale number (n.code+ shade)	Manufacturer	Optical class	Mechanical resistance	Optional requirement	

Euronda Spa © 2019

This catalogue cannot be copied, scanned or reproduced by any means. All rights reserved.

Picture, description, quoting, dimension and technical data have to be considered as purely indicative.

While maintaining the basic characteristics of the product, the manufacturer reserves the right to make modifications at any time with the sole aim of improvement.

Euronda Spa

Via Chizzalunga 1
36066 Sandrigo
Vicenza, Italy
t (+39) 0444 656111
f (+39) 0444 656199
m info@euronda.com
www.euronda.com

Deutschland

Euronda Deutschland GMBH
Am Landwehrbach 5
48341 Altenberge, Deutschland
t (+49) 2505 9389 0
f (+49) 2505 9389 29
m info@euronda.de
www.euronda.de

France

Euronda France Eurl
ZAC les Vallées, Avenue de Bruxelles
60110 Amblainville, France
t + 33 (0) 3 44 06 69 70
f + 33 (0) 3 44 02 03 89
m info@euronda.fr
www.euronda.fr

España

Euronda Spa
Identyd, Av. Cabezo de Torres 4
30007 Murcia, España
t 900102034
m info@euronda.es
www.euronda.es

Russia

OOO Finstar
5/2a Varshavskaya str.
196128 St. Petersburg, Russia
t + 7 (812) 635 88 94
m info@finstartech.ru
www.euronda.com.ru

Euronda®